

Corrugated Tubing

GlobalMed is a leading manufacturer of extruded plastic corrugated tubing. Combined with our injection molded plastic fittings we offer a complete OEM product line of non-invasive disposable and re-usable medical tubing and components for use in anesthesia, respiratory, and other medical and custom applications.

Standard Corrugated Tubing Diameters & Lengths

10mm (I.D.)	15mm (I.D.)	22mm (I.D.)	25mm (I.D.)
6"	4"	4"	40"
12"	6"	6"	48"
24"	12"	12"	60"
36"	24"	24"	72"
40"	40"	40"	96"
48"	60"	48"	120"
60"	72"	60"	192"
72"		72"	
100 & 328 Ft. Coils	100 Foot Coils	100 Foot Coils	--

- Standard 15mm and 22mm medical tubing products are manufactured to ISO specifications
- Coiled tubing is segmented / cuffed in 6" increments or greater
- Standard Tubing Colours: Natural & Translucent Blue
- 6mm, 13mm, 19mm and 32mm profiles are also available

* Custom colours, lengths and profiles are available upon request.

Standard Corrugated Tubing Materials

- EVA (Ethylene Vinyl Acetate)
- PE (Polyethylene)
- PP (Polypropylene)
- Other materials available for custom applications

Our corrugated tubing is used in a variety of medical and custom applications including:

- Anesthesia & Respiratory
- Disposable & Re-usable
- Ventilators
- Sleep Apnea