



Expandable Tubing

GlobalMed is a leading manufacturer of expandable tubing in the Medical Industry with one of the widest size ranges of plastic corrugated and expandable tubing in the market. We offer a complete OEM product line of non-invasive disposable medical tubing and components for use in anesthesia, respiratory, and other medical and custom applications.

Standard Expandable Tubing Diameters & Lengths

10mm (I.D.)	15mm (I.D.)	22mm (I.D.)	38mm (I.D.)
	4"	6"	26"
12"	6"	12"	32"
20"	12"	20"	39"
30"	24"	24"	
40"	40"	30"	
60"	60"	40"	
80"	80"	48"	
120"	108"	60", 70", 80", & 90"	
	120"	108"	
	136"	120"	
		136"	

**Expandable tubing lengths listed are extruded lengths. Due to material design of expandable products there are special length specifications that must be noted. Please contact a member of our Sales Department for more information.

* Custom colours, lengths and profiles are available upon request.

- Standard 15mm and 22mm medical tubing products are manufactured to ISO specifications
- Standard Tubing Colours: Natural & Translucent Blue
- All expandable tubing is packaged in a collapsed state
- Tubing expands to approximately 3x the collapsed length

Expandable Tubing Applications

Our plastic expandable medical tubing is designed for gas flow in a variety of medical and custom applications.

- Anesthesia & Respiratory
- Ventilators
- Disposable & Re-usable