



In addition to our standard corrugated and expandable tubing products GlobalMed also offers 15/22mm corrugated taper tubing in standard lengths ranging from 4 to 72 inches.



15/22mm expandable taper tubing is also available. Although 6 inches is the most common length requested, other lengths are available upon request.



- **Standard Tubing Colours:** Natural & Translucent Blue
- **Standard Materials:** EVA, PE & PP
- **Applications:** - Anesthesia & Respiratory
- Sleep Apnea
- Catheter Mounts

*Custom lengths, colours and material are available upon request.